




BioBrace[®]

Reinforced BioInductive Implant



23 mm x 30 mm

Strengthen Your Repair, Optimize Healing

- Designed for Rotator Cuff Repair and Smaller Tendon Anatomies
- Achieve Native Tendon Strength Within 12 Weeks Through Augmentation^{1,2}
- Provides Up to 515 N of Strength^{1,3}

Features & Benefits

- Provides supplemental strength and load-sharing for up to 2 years^{1,2}
- Generates new tissue via highly porous collagen scaffold^{1,2}
- Fully resorbable^{1,2}
- Robust suture compatibility & handling characteristics designed for open or arthroscopic procedures
- Rapidly absorbs 300% of its weight in blood or other fluids^{1,2}
- Off-the-shelf convenience

¹ K203267, K241906, K242187 - 510(k) Clearance Letter - The BioBrace[®] Implant

² Based on Preclinical Animal Data

³ Based on Internal Testing Data

To learn more, call
1-866-4CONMED (Toll Free)
727-214-3000 (International)
or visit www.CONMED.com.



Imagine what we could do together.




BioBrace® Reinforced BioInductive Implant

Strengthen Your Repair, Optimize Healing

Ordering Information

Description	Catalog Number
BioBrace® Reinforced BioInductive Implant 23mm x 30mm	BB23X30



Call 1-866-4CONMED to schedule a trial with your local CONMED Representative. For more information visit us at www.CONMED.com or scan QR code.

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com