



DuraClip™

Repositionable Hemostasis Clip 11mm & 16mm

For hemostasis within the entire Gastrointestinal (GI) tract, CONMED's DuraClip™ Repositionable Hemostasis Clip is designed to easily and accurately treat defects with unlimited repositionability and improved visualization.




To learn more, call
1-866-4CONMED (Toll Free)
727-214-3310 (International)
or visit www.CONMED.com.



The Next Generation in GI Hemostasis

Smaller Base for Better Visibility

- Up to 60% smaller than competitors



Duraclip™ offers greater compression force than the leading competitor. Duraclip™ is still one of the smallest clips on the market.

Up to 60% shorter post-deployed stem than Boston Res 360.¹

Unlimited Repositioning


- Open and close as needed before deployment.

Optimized Clip Design

- Closer placement and improved visibility.

True 1:1 Rotation

- Precise control in any anatomy.



“We’re getting a very large tissue bite and good apposition with each clip. The short stem doesn’t get in your way, and you can see how well the clips line up and how well that’s closed.”

– Douglas G. Adler, M.D.
GPOEM Case Using DuraClip™ Hemostasis Clips



View Technique




Designed for Defects of Various Sizes

DuraClip™, available in 11mm and 16mm, features a jaw opening that allows physicians to grasp a sizeable amount of tissue – making it a versatile solution.

Treatment indications for DuraClip™ also include:

- Bleeding Ulcers
- Diverticula in the Colon
- As a supplementary method of closure for gastrointestinal tract luminal perforations < 20mm that can be treated conservatively
- Mucosal/sub-mucosal defects less than 3cm and for polyps less than 1.5cm in diameter.



Easy-to-Use Design

The CONMED handle has a familiar, intuitive design for one or two-hand clip deployment that gives both a tactile and audible response once clip is deployed.




Precise Positioning

- Achieve accurate and precise rotation of clip
- Flexible distal end enables rotation over ERCP scope elevator
- Accommodates insertion and rotation in retroflexed scope

Easy Deployment

- Design enables simple one or two handle clip deployment
- Unique, small post-deployed clip base improves field of view
- Tactile and audible response confirm deployment

Duraclip is available in 11mm and 16mm jaw width.



The 11mm Duraclip™ has a 60% shorter clip base/stem than Res 360,² and the 16mm Duraclip™ has a 39% shorter clip base/stem.





DuraClip™

Repositionable Hemostasis Clip 11mm & 16mm

Ordering Information

Description	Sheath Length	Min Channel Diameter	Jaw Opening	Qty / Case	Catalog Number
DuraClip™ Repositionable Hemostasis Clip	235cm	2.8mm	11mm	10 per box	DC0235
DuraClip™ Repositionable Hemostasis Clip	235cm	2.8mm	16mm	10 per box	DC0235W

*Data on file: TR16-756 and TR17-0300. Comparison is to Boston Scientific Resolution™ Clip.

**Data on file: TR16-757 and TR17-0300.

¹ When comparing DuraClip™ to Boston Scientific Revolution 360. Data on file, PDD2324060 .





DuraClip™ Repositionable Hemostasis Clip



CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | www.CONMED.com



Service & Support