



PlumeSafe[®] X5[™]

SMOKE MANAGEMENT SYSTEM


PlumeSafe[®] X5[™] offers enhanced efficiency and safety through its streamlined, user-friendly design and advanced 5-stage filter—the first of its kind. Welcome to the next generation of smoke management systems.

FIRST OF ITS KIND 5-STAGE FILTER




USER FRIENDLY

Easily select between 3 procedure modes with a simple touch-sensitive display.



ADVANCED FILTRATION

5-stage filter that captures particles down to .01 microns at 99.9997% efficiency.*



OPTIMIZED DESIGN

Dual 7/8" ports to connect two smoke pencils without any additional accessories.

EXCEPTIONAL SMOKE EVACUATION **STARTS HERE**

* Data on file.



FACILITATE A HEALTHIER ENVIRONMENT FOR **EVERYONE**
IN YOUR OPERATING ROOM WITH PLUMESAFE® X5™.

Connect with your local representative today or call **1-866-4CONMED**.

ORDERING & PRODUCT INFORMATION

DESCRIPTION	CATALOG NUMBER
CONMED PlumeSafe® X5™ Smoke Management System, 120V	PX5-120
CONMED PlumeSafe® X5™ Smoke Management System, 220V	PX5-220
CONMED PlumeSafe® X5™ Smoke Management System, 220V, Sealed	PX5-220P1
CONMED PlumeSafe® X5™ Filter, 1/Case, 22mm	PX5-135-22
CONMED PlumeSafe® X5™ Filter, 4/Case, 22mm	PX5-1354-22



CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com